

## Anti-Alkaline phosphatase/ALPL Antibody Picoband®

Catalog Number: A01008

### About ALPL

Alkaline phosphatase (ALPL) removes phosphate groups from the 5' end of DNA and RNA, and from proteins, at high pH. Most mammals have 4 different isozymes: placental, placental like, intestinal and non tissue specific (found in liver, kidney and bone). Tissues with particularly high concentrations of ALP include the liver, bile ducts, placenta, and bone. ALPL is the alkaline phosphatase of skin fibroblasts, the tissue-nonspecific type, and that it is active toward millimolar concentrations of the putative natural substrates phosphoethanolamine (PEA) and pyridoxal-5-prime-phosphate (PLP). ALPL gene exists in single copy in the haploid genome and is composed of 12 exons distributed over more than 50 kb. Damaged or diseased tissue releases enzymes into the blood, so serum ALP measurements can be abnormal in many conditions, including bone disease and liver disease.

### Overview


Product Name	Anti-Alkaline phosphatase/ALPL Antibody Picoband®
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Alkaline phosphatase/ALPL Antibody Picoband® catalog # A01008. Tested in ELISA, IF, ICC, WB applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.01mg NaN <sub>3</sub> .
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05186

### Technical Details


Immunogen	E.coli-derived human Alkaline phosphatase/ALPL recombinant protein (Position: K229-H482).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot, and HRP Conjugated anti-Rabbit IgG Super Vision Assay Kit (SV0002-1) for ICC.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG

Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5ug/ml, Human, Mouse, Rat Immunocytochemistry/Immunofluorescence, 5ug/ml, Human ELISA, 0.1-0.5ug/ml, -

## Anti-Alkaline phosphatase/ALPL Antibody Picoband® (A01008) Images



IF analysis of ALPL using anti-ALPL antibody (A01008). ALPL was detected in immunocytochemical section of A431 cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5ug/mL rabbit anti-ALPL Antibody (A01008) overnight at 4°C. DyLight®594 Conjugated Goat Anti-Rabbit IgG (BA1142) was used as secondary antibody at 1:100 dilution and incubated for 30 minutes at 37°C. The section was counterstained with DAPI. Visualize using a fluorescence microscope and filter sets appropriate for the label used.



Western blot analysis of ALPL using anti-ALPL antibody (A01008). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 50ug of sample under reducing conditions. Lane 1: human HeLa whole cell lysates, Lane 2: human HEK293 whole cell lysates, Lane 3: rat kidney tissue lysates, Lane 4: rat RH-35 whole cell lysates, Lane 5: mouse kidney tissue lysates, Lane 6: mouse HEPA1-6 whole cell lysates. After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes. Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-ALPL antigen affinity purified polyclonal antibody (Catalog # A01008) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for ALPL at approximately 80KD. The expected band size for ALPL is at 57KD.

## 2 Publications Citing This Product

1. PubMed ID: 32018039, Chen G,Wang Q,Li Z,Yang Q,Liu Y,Du Z,Zhang G,Song Y.Circular RNA CDR1as promotes adipogenic and suppresses osteogenic differentiation of BMSCs in steroid-induced osteonecrosis of the femoral head.Bone.2020 Apr;133:115258.doi:10.1016/j.bone.2020.115258.Epub

2. PubMed ID: 32857019, Zhang Y,Yin C,Cheng Y,Huang X,Liu K,Cheng G,Li Z.Electrospinning Nanofiber-Reinforced Aerogels for the Treatment of Bone Defects.Adv Wound Care (New Rochelle).2020 Aug;9(8):441-452.doi:10.1089/wound.2018.0879.Epub 2019 Jun 5.PMID:32857019;PMCID:PMC7382394

Visit [bosterbio.com/anti-alkaline-phosphatase-alpl-picoband-trade-antibody-a01008-boster.html](http://bosterbio.com/anti-alkaline-phosphatase-alpl-picoband-trade-antibody-a01008-boster.html) to see all 2 publications.

Submit a product review to [Biocompare.com](http://Biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Alkaline phosphatase/ALPL Antibody

For Research Use Only. Not for use in diagnostic procedures.