

Anti-Integrin Beta 4 ITGB4 Antibody

Catalog Number: A01015


Overview

Product Name	Anti-Integrin Beta 4 ITGB4 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Integrin Beta 4 ITGB4 Antibody catalog # A01015. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16144


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 IHC: 1:50-1:200

Anti-Integrin Beta 4 ITGB4 Antibody (A01015) Images



Western blot analysis of Integrin beta 4 on Jurkat cells lysates using anti-Integrin beta 4 antibody at 1/500 dilution. Positive control □ Line 1: A431 Line 2: A549



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue using anti-Integrin beta 4 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin Beta 4 ITGB4 Antibody

For Research Use Only. Not for use in diagnostic procedures.