

Anti-Glutathione peroxidase 1 GPX1 Antibody

Catalog Number: A01019-1

Overview

Product Name	Anti-Glutathione peroxidase 1 GPX1 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Glutathione peroxidase 1 GPX1 Antibody catalog # A01019-1. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07203

Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-Glutathione peroxidase 1 GPX1 Antibody (A01019-1) Images




Figure 2. Immunocytochemistry staining of GPX1 using Anti-Glutathione peroxidase 1 GPX1 Antibody (A01019-1).

ICC staining Glutathione Peroxidase 1 in 293 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 1. Western blotting validation for Anti-Glutathione peroxidase 1 GPX1 Antibody A01019-1

Western blot analysis of Glutathione Peroxidase 1 on THP-1 cells lysates using anti-Glutathione Peroxidase 1 antibody at 1/1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glutathione peroxidase 1 GPX1 Antibody