

Anti-FOXA2 Antibody

Catalog Number: A01032-1


Overview

Product Name	Anti-FOXA2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-FOXA2 Antibody catalog # A01032-1. Tested in WB,ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y261

Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200

Anti-FOXA2 Antibody (A01032-1) Images



Western blot analysis of FOXA2 on HT-29 cells lysates using anti-FOXA2 antibody at 1/500 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FOXA2 Antibody

For Research Use Only. Not for use in diagnostic procedures.