

Anti-NeuroD1 Antibody

Catalog Number: A01038


Overview

Product Name	Anti-NeuroD1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NeuroD1 Antibody catalog # A01038. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13562

Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IP: 1:10-1:50

Anti-NeuroD1 Antibody (A01038) Images



Western blot analysis of NeuroD1 on different cells lysates using anti-NeuroD1 antibody at 1/500 dilution. Positive control: Line 1: human brain Line 2:SH-SY5Y

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NeuroD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.