

Anti-NEUROD1 Antibody Picoband® PE Conjugated

Catalog Number: A01038-2-PE

About NEUROD1

This gene encodes a member of the NeuroD family of basic helix-loop-helix (bHLH) transcription factors. The protein forms heterodimers with other bHLH proteins and activates transcription of genes that contain a specific DNA sequence known as the E-box. It regulates expression of the insulin gene, and mutations in this gene result in type II diabetes mellitus.

Overview

| | |
|----------------------|--|
| Product Name | Anti-NEUROD1 Antibody Picoband® PE Conjugated |
| Reactive Species | Human, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q13562 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the N-terminus of human NEUROD1. Human NEUROD1 shares 90.9% amino acid (aa) sequence identity with both mouse and rat NEUROD1. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NEUROD1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.