

## Anti-CD79a Antibody

Catalog Number: A01047-2


### Overview

Product Name	Anti-CD79a Antibody
Reactive Species	Human
Description	Boster Bio Anti-CD79a Antibody catalog # A01047-2. Tested in WB,IF,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11912


### Technical Details

Immunogen	Recombinant protein within human CD79a aa 33-143.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified.
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 IF: 1:50-1:100 FC: 1:50-1:100

## Anti-CD79a Antibody (A01047-2) Images



Western blot analysis of CD79a on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature. Positive control: Lane 1: Daudi cell lysate Lane 2: Raji cell lysate



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-CD79a antibody. The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 8.0-8.4) for 20 minutes. The tissues were blocked in 5% BSA for 30 minutes at room temperature, washed with ddH<sub>2</sub>O and PBS, and then probed with ET7109-33 at 1/100 dilution, for 30 minutes at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. Counter stained with hematoxylin and mounted with DPX.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-CD79a Antibody

For Research Use Only. Not for use in diagnostic procedures.