

Anti-Phospho-LATS1/2 (Thr1079/1041) Antibody

Catalog Number: A01051T1079-1


Overview

Product Name	Anti-Phospho-LATS1/2 (Thr1079/1041) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-LATS1/2 (Thr1079/1041) Antibody catalog # A01051T1079-1. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95835

Technical Details

Immunogen	Synthesized phospho derived from human LATS1/2 (Phospho-Thr1079/1041)
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:10000-20000

Anti-Phospho-LATS1/2 (Thr1079/1041) Antibody (A01051T1079-1) Images



Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-LATS1/2 (Thr1079/1041) Antibody

For Research Use Only. Not for use in diagnostic procedures.