

Anti-YB1/Ybx1 Antibody Picoband® Biotin Conjugated

Catalog Number: A01054-1-Biotin

About Ybx1

Y box binding protein 1 also known as Y-box transcription factor or nuclease-sensitive element-binding protein 1 is a protein that in humans is encoded by the YBX1 gene. Enables chromatin binding activity; mRNA binding activity; and single-stranded DNA binding activity. Involved in epidermis development and negative regulation of cellular senescence. Acts upstream of or within several processes, including cellular response to interleukin-7; negative regulation of striated muscle cell differentiation; and negative regulation of transcription by RNA polymerase II. Located in cytoplasm and nucleus. Part of histone pre-mRNA 3'end processing complex. Is expressed in several structures, including branchial arch; central nervous system; facial prominence; gonad; and retina. Orthologous to human YBX1 (Y-box binding protein 1).

Overview

Product Name	Anti-YB1/Ybx1 Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P62960

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse YB1/Ybx1, identical to the related rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users.

Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.
ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-YB1/Ybx1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.