

Anti-RSK1 p90 (Phospho-S380) RPS6KA1 Antibody

Catalog Number: A01058S380-2


Overview

Product Name	Anti-RSK1 p90 (Phospho-S380) RPS6KA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RSK1 p90 (Phospho-S380) RPS6KA1 Antibody catalog # A01058S380-2. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15418


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser380 of human RSK1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200

Anti-RSK1 p90 (Phospho-S380) RPS6KA1 Antibody (A01058S380-2) Images



Western blot analysis of Phospho-RSK1(S380) on A431 cell lysates using anti-Phospho-RSK1(S380) antibody at 1/1,000 dilution.



ICC staining Phospho-RSK1(S380) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RSK1 p90 (Phospho-S380) RPS6KA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.