

Anti-Rsk-1 (phospho-T359/S363) RPS6KA1 Antibody

Catalog Number: A01058T359S363


Overview

Product Name	Anti-Rsk-1 (phospho-T359/S363) RPS6KA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Rsk-1 (phospho-T359/S363) RPS6KA1 Antibody catalog # A01058T359S363. Tested in WB,IP,IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15418


Technical Details

Immunogen	Synthetic phosphopeptide derived from human Rsk-1 around the phosphorylation site of Thr359/Ser363.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200 IP: 1:50-1:200

Anti-Rsk-1 (phospho-T359/S363) RPS6KA1 Antibody (A01058T359S363) Images




Western blot (WB) analysis of p-Rsk-1 (T359/S363) pAb at 1:2000 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:HEK293T treated with EGF(100ng/ml,30 minutes) whole cell lysate(40ug)




Immunofluorescence analysis of CT26 cells using Rsk-1 (phospho-T359/S363) antibody at dilution of 1:50.

WB:



Immunoprecipitation of HepG2 cell lysate using p-Rsk-1 (T359/S363) pAb (Sepharose Bead Conjugate) #BD0048(lane 2 and lane 3) and Nonspecific IgG Control (Sepharose Bead Conjugate) #BD0048 (lane 4) .Lane 1 is 30% input. The western blot was probed using p-Rsk-1 (T359/S363) #BS4276. “↑”[supernatant]; “↓”(deposition)



Immunohistochemistry of paraffin-embedded Rat Brain using Rsk-1 (phospho-T359/S363) antibody at dilution of 1:50.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rsk-1 (phospho-T359/S363) RPS6KA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.