

## Anti-RSK1 p90 (Phospho-T359/S363) RPS6KA1 Antibody

Catalog Number: A01058T359S363-1


### Overview

Product Name	Anti-RSK1 p90 (Phospho-T359/S363) RPS6KA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RSK1 p90 (Phospho-T359/S363) RPS6KA1 Antibody catalog # A01058T359S363-1. Tested in WB,ICC,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15418


### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr359 and Ser363 of human RSK1 p90
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

## Anti-RSK1 p90 (Phospho-T359/S363) RPS6KA1 Antibody (A01058T359S363-1) Images



Western blot analysis of Phospho-RSK1 p90(T359+S363) on different lysates using anti-Phospho-RSK1 p90(T359+S363) antibody at 1/1,000 dilution. Positive control: Lane 1: A431  
Lane 2: 293 Lane 3: MCF-7



ICC staining Phospho-RSK1 p90(T359+S363) in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RSK1 p90 (Phospho-T359/S363) RPS6KA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.