

Anti-Blood group Rh(D) polypeptide RHD Antibody

Catalog Number: A01082-1


Overview

Product Name	Anti-Blood group Rh(D) polypeptide RHD Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Blood group Rh (D) polypeptide RHD Antibody catalog # A01082-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q02161

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 168-213 of Human RHD.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Blood group Rh(D) polypeptide RHD Antibody (A01082-1) Images



Western blot (WB) analysis of RHD polyclonal antibody at 1:500 dilution Lane1:The Brain tissue lysate of Rat(30ug) Lane2:The Brain tissue lysate of Mouse(30ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:H9C2 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Blood group Rh(D) polypeptide RHD Antibody

For Research Use Only. Not for use in diagnostic procedures.