

Anti-Stat5a/b (A688) Antibody

Catalog Number: A01087


Overview

Product Name	Anti-Stat5a/b (A688) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Stat5a/b (A688) Antibody catalog # A01087. Tested in WB,IF applications. This antibody reacts with Human,Mouse.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42229/P51692


Technical Details

Immunogen	Recombinant protein.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IF: 1:50-1:200

Anti-Stat5a/b (A688) Antibody (A01087) Images



Western blot (WB) analysis of Stat5a/b pAb at 1:500 dilution
Lane1:PC3 whole cell lysate(40ug) Lane2:CT-26 whole cell lysate(40ug) Lane3:THP-1 whole cell lysate(40ug)
Lane4:MCF-7 whole cell lysate(40ug)



Immunofluorescence analysis of MCF-7 cells using Stat5a/b (A688) antibody at dilution of 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat5a/b (A688) Antibody

For Research Use Only. Not for use in diagnostic procedures.