

Anti-Stat5a/b (phospho-Y694) Antibody

Catalog Number: A01087Y694


Overview

Product Name	Anti-Stat5a/b (phospho-Y694) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Stat5a/b (phospho-Y694) Antibody catalog # A01087Y694. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42229/P51692

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Stat5a around the phosphorylation site of Tyrosine 694.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-Stat5a/b (phospho-Y694) Antibody (A01087Y694) Images



Immunohistochemistry (IHC) analyzes of p-Stat5a/b (Y694) pAb in paraffin-embedded human breast carcinoma tissue

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat5a/b (phospho-Y694) Antibody

For Research Use Only. Not for use in diagnostic procedures.