

## Anti-Stat5 (Phospho-Y694) STAT5A Antibody

Catalog Number: A01087Y694-1

### Overview


Product Name	Anti-Stat5 (Phospho-Y694) STAT5A Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Stat5 (Phospho-Y694) STAT5A Antibody catalog # A01087Y694-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42229

### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Tyr694 of human Stat5.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 IHC: 1:50-1:200

## Anti-Stat5 (Phospho-Y694) STAT5A Antibody (A01087Y694-1) Images

---



Western blot analysis of Phospho-Stat5(Y694) on mouse liver lysates using anti-Phospho-Stat5(Y694) antibody at 1/1,000 dilution.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat5 (Phospho-Y694) STAT5A Antibody

For Research Use Only. Not for use in diagnostic procedures.