

Anti-C4a/b Antibody

Catalog Number: A01095


Overview

Product Name	Anti-C4a/b Antibody
Reactive Species	Human
Description	Boster Bio Anti-C4a/b Antibody catalog # A01095. Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	POCOL4

Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human C4a/b.
Isotype	IgG
Form	Liquid
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB, 1:500-1:2000 IHC-P, 1:100-1:300</p> <p>For protocols please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>


Anti-C4a/b Antibody (A01095) Images



Western Blot (WB) analysis of AD293 cells using C4a/b Polyclonal antibody.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Liver, antibody was diluted at 1:100.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain, antibody was diluted at 1:100.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Liver, antibody was diluted at 1:100.

1 Publications Citing This Product

1. PubMed ID: 2617767, Alternative Complement Activity in the Egg Cytosol of Amphioxus Branchiostoma belcheri: Evidence for the Defense Role of Maternal Complement Components

Visit bosterbio.com/anti-c4a-b-antibody-a01095-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C4a/b Antibody