

Anti-C4a/b Antibody

Catalog Number: A01095


Overview

Product Name	Anti-C4a/b Antibody
Reactive Species	Human
Description	Boster Bio Anti-C4a/b Antibody catalog # A01095. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	POC0L4

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human C4A/C4B. AA range:721-770
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:10000 IF 1:50-200


Anti-C4a/b Antibody (A01095) Images



Western Blot (WB) analysis of AD293 cells using C4a/b Polyclonal antibody.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Liver, antibody was diluted at 1:100.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain, antibody was diluted at 1:100.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Liver, antibody was diluted at 1:100.

1 Publications Citing This Product

1. PubMed ID: 2617767, Alternative Complement Activity in the Egg Cytosol of Amphioxus Branchiostoma belcheri: Evidence for the Defense Role of Maternal Complement Components

Visit bosterbio.com/anti-c4a-b-antibody-a01095-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C4a/b Antibody

For Research Use Only. Not for use in diagnostic procedures.