

## Anti-C4 Antibody

Catalog Number: A01095-2


### Overview

Product Name	Anti-C4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-C4 Antibody catalog # A01095-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P0COL4


### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:200-1:500 ICC: 1:50-1:200 IHC: 1:50-1:200

## Anti-C4 Antibody (A01095-2) Images



Western blot analysis of C4 on human liver tissue lysate using anti-C4 antibody at 1/1,000 dilution.



ICC staining C4 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-C4 Antibody

For Research Use Only. Not for use in diagnostic procedures.