

Anti-C4 Antibody

Catalog Number: A01095-2

Overview

Product Name	Anti-C4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-C4 Antibody catalog # A01095-2. Tested in IF, IHC, ICC, WB applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P0COL4

Technical Details

Form	Liquid
Purification	ProA affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB:1:200-1:500 ICC:1:50-1:200 IHC:1:50-1:200</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-C4 Antibody (A01095-2) Images




Figure 2. Immunocytochemistry staining of C4A using Anti-C4 Antibody (A01095-2).

ICC staining C4 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 1. Western blotting validation for Anti-C4 Antibody A01095-2

Western blot analysis of C4 on human liver tissue lysate using anti-C4 antibody at 1/1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C4 Antibody