

Anti-IL-23 Antibody

Catalog Number: A01097

About IL23A

Like interleukin-27 (IL-27), IL-23 is a recently discovered member of the IL-6/IL-12 family of proinflammatory and immunoregulatory cytokines. It exists as a heterodimer composed of the IL-12p40 subunit and a novel p19 subunit. IL-23 is secreted by activated dendritic cells, macrophages, and monocytes. Its biological activities include enhancing the proliferation of memory T cells and the production of IFN-gamma, IL-12, and TNF-alpha from activated T cells, and can stimulate macrophages to produce TNF-alpha and nitric oxide. It has also been shown to possess potent anti-tumor and anti-metastatic activity in mouse models of cancer, suggesting a potential role for IL-23 in therapeutic treatment of cancer.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-IL-23 Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL-23 Antibody (Catalog # A01097). Tested in ELISA, WB, IHC-P, ICC, IF applications. This antibody reacts with Human.
Application	ELISA, IF, IHC-P, ICC, WB
Clonality	Polyclonal
Formulation	IL-23 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	IL-23 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q9NPF7

Technical Details


Immunogen	Anti-IL-23 antibody was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of mature human IL-23. The immunogen is located the last 50 amino acids of IL-23.
Predicted Reactive Species	Guinea Pig, Horse, Mouse, Rat
Isotype	IgG

Form	Liquid
Concentration	1 mg/mL
Purification	IL-23 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	WB: 1-4 ug/mL; ICC: 10 ug/mL; IF: 20 ug/mL, IHC: 5 ug/mL. Antibody validated: Western Blot, Immunocytochemistry, Immunofluorescence and Immunohistochemistry in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-IL-23 Antibody (A01097) Images




Western Blot Validation in Human Raji Cells Loading: 15 ug of lysates per lane. Antibodies: IL-23 A01097 (1 ug/mL), 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution. Lane A: Absence of blocking peptide Lane B: Presence of blocking peptide.



Western Blot Validation in Human Cell Lines Loading: 20 ug of lysates per lane. Antibodies: IL-23 A01097, (2 ug/mL), 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Western Blot Validation with Recombinant Protein Loading: 30 ng of human IL-23 recombinant protein. Antibodies: IL-23 A01097 (4 ug/mL), 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Immunofluorescence Validation of IL-23 in Raji Cells
Immunofluorescent analysis of 4% paraformaldehyde-fixed Raji Cells labeling IL-23 with A01097 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (red).


Immunocytochemistry Validation of IL-23 in Raji Cells
Immunocytochemical analysis of Raji cells using anti-IL-23 antibody (A01097) at 10 ug/ml. Cells were fixed with formaldehyde and blocked with 10% serum for 1 h at RT; antigen retrieval was by heat mediation with a citrate buffer (pH6). Samples were incubated with primary antibody overnight at 4° C. A goat anti-rabbit IgG H&L (HRP) at 1/250




was used as secondary. Counter stained with Hematoxylin.



Immunofluorescence Validation of IL-23 in Human Pancreas Tissue
Immunofluorescent analysis of 4% paraformaldehyde-fixed human pancreas tissue labeling IL-23 with A01097 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (red) and DAPI staining (blue).



Immunohistochemistry Validation of IL-23 in Human Pancreas Tissue
Immunohistochemical analysis of paraffin-embedded human pancreas tissue using anti-IL-23 antibody (A01097) at 5 ug/ml. Tissue was fixed with formaldehyde and blocked with 10% serum for 1 h at RT; antigen retrieval was by heat mediation with a citrate buffer (pH6). Samples were incubated with primary antibody overnight at 4° C. A goat anti-rabbit IgG H&L (HRP) at 1/250 was used as secondary. Counter stained with Hematoxylin.



Immunohistochemistry Validation of IL-23 in Human Abdominal Aortic Aneurysm (AAA) Tissues (Yan et al., 2019)

Anti-IL-23 Antibody

For Research Use Only. Not for use in diagnostic procedures.