

Anti-Liver Arginase ARG1 Antibody

Catalog Number: A01106-2


Overview

Product Name	Anti-Liver Arginase ARG1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Liver Arginase ARG1 Antibody catalog # A01106-2. Tested in WB,IHC,IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05089


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 IP: 1:10-1:50

Anti-Liver Arginase ARG1 Antibody (A01106-2) Images



Western blot analysis of Liver Arginase on human liver tissue lysate using anti-Liver Arginase antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-Liver Arginase antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Liver Arginase ARG1 Antibody

For Research Use Only. Not for use in diagnostic procedures.