

Anti-Leucyl tRNA synthetase /LARS Antibody Picoband® Biotin Conjugated

Catalog Number: A01109-Biotin

About LARS

Leucyl-tRNA synthetase, cytoplasmic is an enzyme that in humans is encoded by the LARS gene. This gene is mapped to 5q32. It encodes a cytosolic leucine-tRNA synthetase, a member of the class I aminoacyl-tRNA synthetase family. The encoded enzyme catalyzes the ATP-dependent ligation of L-leucine to tRNA (Leu). It is found in the cytoplasm as part of a multisynthetase complex and interacts with the arginine tRNA synthetase through its C-terminal domain. A mutation in this gene was found in affected individuals with infantile liver failure syndrome 1. Alternatively spliced transcript variants of this gene have been observed.

Overview

Product Name	Anti-Leucyl tRNA synthetase /LARS Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9P2J5

Technical Details

Immunogen	E. coli-derived human Leucyl tRNA synthetase/LARS recombinant protein (Q628-Q855).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Leucyl tRNA synthetase /LARS Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.