

## Anti-Dopamine beta Hydroxylase/DBH Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A01110-Fluoro647

### About DBH

Dopamine beta-hydroxylase (DBH), also known as dopamine beta-monoxygenase, is an enzyme (EC1.14.17.1) that in humans is encoded by the DBH gene. Dopamine beta-hydroxylase catalyzes the chemical reaction. It is mapped to 9q34.2. The protein encoded by this gene is an oxidoreductase belonging to the copper type II, ascorbate-dependent monoxygenase family. The encoded protein, expressed in neurosecretory vesicles and chromaffin granules of the adrenal medulla, catalyzes the conversion of dopamine to norepinephrine, which functions as both a hormone and as the main neurotransmitter of the sympathetic nervous system. The enzyme encoded by this gene exists in both soluble and membrane-bound forms, depending on the absence or presence, respectively, of a signal peptide. Mutations in this gene cause dopamine beta-hydroxylase deficiency in human patients, characterized by deficits in autonomic and cardiovascular function, including hypotension and ptosis. Polymorphisms in this gene may play a role in a variety of psychiatric disorders.

### Overview

Product Name	Anti-Dopamine beta Hydroxylase/DBH Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P09172

### Technical Details

Immunogen	E.coli-derived human Dopamine beta Hydroxylase/DBH recombinant protein (Position: S40-N545).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.

Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 30315688, Ren S,Chen M,Yang L,Liu Z.5-Hydroxytryptamine and Dopamine Neurons in the Cerebellum of the New-Hatching Yangtze Alligator Alligator sinensis.Anat Rec(Hoboken).2019 Jun;302(6):861-868.doi:10.1002/ar.23982.Epub 2018 Nov 13.PMID:30315688.

Visit [bosterbio.com/anti-dopamine-beta-hydroxylase-dbh-picoband-trade-antibody-a01110-boster.html](https://bosterbio.com/anti-dopamine-beta-hydroxylase-dbh-picoband-trade-antibody-a01110-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dopamine beta Hydroxylase/DBH Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.