

## Anti-Podoplanin PDPN Antibody

Catalog Number: A01124-1


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-Podoplanin PDPN Antibody  |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-Podoplanin PDPN Antibody catalog # A01124-1. Tested in WB,IHC applications. This antibody reacts with Human.               |
| Application          | IHC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q86YL7   |


### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | recombinant protein                      |
| Cross Reactivity    | No cross reactivity with other proteins. |
| Isotype             | IgG                                      |
| Form                | Liquid                                   |
| Concentration       | 1 mg/ml                                  |
| Purification        | ProA affinity purified                   |
| Suggested Dilutions | WB: 1:2,000-1:5,000<br>IHC: 1:50-1:200   |

## Anti-Podoplanin PDPN Antibody (A01124-1) Images



Western blot analysis of Podoplanin on human lung cells lysates using anti-Podoplanin antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human lung tissue using anti-Podoplanin antibody. Counter stained with hematoxylin.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Podoplanin PDPN Antibody

For Research Use Only. Not for use in diagnostic procedures.