

Anti-IKB Alpha NFKBIA Antibody

Catalog Number: A01139


Overview

Product Name	Anti-IKB Alpha NFKBIA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IKB Alpha NFKBIA Antibody catalog # A01139. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P25963


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-IKB Alpha NFKBIA Antibody (A01139) Images



Western blot analysis of IKB alpha on different lysates using anti-IKB alpha antibody at 1/1,000 dilution. Positive control:
Lane 1: PC-12 Lane 2: NIH/3T3 Lane 3: Hela



ICC staining IKB alpha in Hela cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 10.1007/s00253-021-11719-8, Surfactin-oleogel with therapeutic potential for inflammatory acne vulgaris induced by *Propionibacterium acnes*

Visit bosterbio.com/anti-ikb-alpha-nfkb-antibody-a01139-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IKB Alpha NFKBIA Antibody

For Research Use Only. Not for use in diagnostic procedures.