

Anti-IL-1 Alpha (F8) IL1A Antibody

Catalog Number: A01144-2

About IL1A

This protein mediates the voltage-dependent sodium ion permeability of excitable membranes. Assuming opened or closed conformations in response to the voltage difference across the membrane, the protein forms a sodium-selective channel through which Na⁺ ions may pass in accordance with their electrochemical gradient. This sodium channel may be present in both denervated and innervated skeletal muscle.

George A.L. Jr., Ann. Neurol. 31:131-137(1992).

Wang J., Biochem. Biophys. Res. Commun. 182:794-801(1992).

Tsujino A., Proc. Natl. Acad. Sci. U.S.A. 100:7377-7382(2003).

Overview


Product Name	Anti-IL-1 Alpha (F8) IL1A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IL-1 Alpha (F8) IL1A Antibody catalog # A01144-2. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01583

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 91-140 of Human IL-1alpha.
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml

Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>IHC: 1:50-1:200</p>

Anti-IL-1 Alpha (F8) IL1A Antibody (A01144-2) Images



Immunohistochemistry (IHC) analyzes of IL-1alpha (F8) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-1 Alpha (F8) IL1A Antibody