

## Anti-CD27 antigen CD27 Antibody

Catalog Number: A01148-1


### Overview

Product Name	Anti-CD27 antigen CD27 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-CD27 antigen CD27 Antibody catalog # A01148-1. Tested in WB,ICC,IHC,Flow Cytometry applications. This antibody reacts with Human,rat.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P26842


### Technical Details

Immunogen	Intracellular domain corresponding to Human CD27.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 ICC: 1:50-1:200 FC: 1:50-1:100

## Anti-CD27 antigen CD27 Antibody (A01148-1) Images



Western blot analysis of CD27 on Raji cell using anti-CD27 antibody at 1/500 dilution.



ICC staining CD27 in K562 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD27 antigen CD27 Antibody

For Research Use Only. Not for use in diagnostic procedures.