

Anti-PBR TSPO Antibody

Catalog Number: A01153-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-PBR TSPO Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-PBR TSPO Antibody catalog # A01153-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P30536 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 124-169 of Human PBR. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000 |

Anti-PBR TSPO Antibody (A01153-1) Images



Western blot (WB) analysis of PBR pAb at 1:2000 dilution
Lane1:PMVEC whole cell lysate(40ug) Lane2:HepG2 whole cell lysate(40ug) Lane3:3T3-L1 whole cell lysate(40ug)
Lane4:A2780 whole cell lysate(40ug) Lane5:A549 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PBR TSPO Antibody

For Research Use Only. Not for use in diagnostic procedures.