

Anti-IL-11 Monoclonal Antibody

Catalog Number: A01160


Overview

Product Name	Anti-IL-11 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL-11 Monoclonal Antibody catalog # A01160. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 2D3
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P20809

Technical Details

Immunogen	Purified recombinant fragment of human IL-11 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-IL-11 Monoclonal Antibody (A01160) Images



Western Blot analysis using IL-11 Monoclonal Antibody against truncated IL-11 recombinant protein.

1 Publications Citing This Product

1. PubMed ID: 10.1007/s12195-019-00606-y, Cas9 Protein Triggers Differential Expression of Inherent Genes Especially NGFR Expression in 293T Cells

Visit bosterbio.com/anti-il-11-antibody-a01160-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-11 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.