

Anti-PLGF PGF Antibody

Catalog Number: A01164-3


Overview

Product Name	Anti-PLGF PGF Antibody
Reactive Species	Human
Description	Boster Bio Anti-PLGF PGF Antibody catalog # A01164-3. Tested in WB,ICC,IHC applications. This antibody reacts with Human.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49763


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-PLGF PGF Antibody (A01164-3) Images



Western blot analysis of PLGF on SH-SY5Y cell lysate using anti-PLGF antibody at 1/1,000 dilution.



ICC staining PLGF in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLGF PGF Antibody

For Research Use Only. Not for use in diagnostic procedures.