

## Anti-CD201 PROCR Antibody

Catalog Number: A01166


### Overview

Product Name	Anti-CD201 PROCR Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD201 PROCR Antibody catalog # A01166. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UNN8


### Technical Details

Immunogen	Synthesized peptide derived from Endothelial protein C receptor at AA range: 141-190
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-300 ELISA 1:20000 IF 1:50-200


## Anti-CD201 PROCR Antibody (A01166) Images



Immunohistochemical analysis of paraffin-embedded mouse spleen, antibody was diluted at 1:100



Western blot analysis of Rat-barin lysis using CD201 antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human spleen, antibody was diluted at 1:100

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD201 PROCR Antibody

For Research Use Only. Not for use in diagnostic procedures.