

Anti-Transcriptional activator GLI3 Gli3 Antibody

Catalog Number: A01171-1


Overview

| | |
|----------------------|---|
| Product Name | Anti-Transcriptional activator GLI3 Gli3 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-Transcriptional activator GLI3 Gli3 Antibody catalog # A01171-1. Tested in WB, ICC applications. This antibody reacts with Human. |
| Application | ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P10071 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant protein within human Gli3 aa 1300-1600. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:500-1:2,000 ICC: 1:50-1:100 |

Anti-Transcriptional activator GLI3 Gli3 Antibody (A01171-1) Images



Western blot analysis of Gli3 on PC-3M cell lysate using anti-Gli3 antibody at 1/1,000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcriptional activator GLI3 Gli3 Antibody

For Research Use Only. Not for use in diagnostic procedures.