

Anti-Cytochrome P450 2C8 CYP2C8 Antibody

Catalog Number: A01179

Overview

Product Name	Anti-Cytochrome P450 2C8 CYP2C8 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cytochrome P450 2C8 CYP2C8 Antibody catalog # A01179. Tested in ELISA, IF, IHC-P applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10632

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human CYP2C8. AA range:241-290
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC: 1:100-1:300 IF 1:50-200 ELISA 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 2C8 CYP2C8 Antibody

For Research Use Only. Not for use in diagnostic procedures.