

Anti-Cytochrome P450 2C8 CYP2C8 Rabbit Polyclonal Antibody

Catalog Number: A01179-1

Overview

Product Name	Anti-Cytochrome P450 2C8 CYP2C8 Rabbit Polyclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cytochrome P450 2C8 CYP2C8 Rabbit Polyclonal Antibody catalog # A01179-1. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10632

Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human CYP2C8.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 2C8 CYP2C8 Rabbit Polyclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.