

## Anti-Choline Acetyltransferase/Chat Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A01192-5-Fluoro488

### About Chat

Choline acetyltransferase (commonly abbreviated as ChAT, but sometimes CAT) is a transferase enzyme responsible for the synthesis of the neurotransmitter acetylcholine. In humans, the choline acetyltransferase enzyme is encoded by the CHAT gene. This gene product is a characteristic feature of cholinergic neurons, and changes in these neurons may explain some of the symptoms of Alzheimer's disease. Polymorphisms in this gene have been associated with Alzheimer's disease and mild cognitive impairment. Mutations in this gene are associated with congenital myasthenic syndrome associated with episodic apnea. Multiple transcript variants encoding different isoforms have been found for this gene, and some of these variants have been shown to encode more than one isoform.

### Overview

Product Name	Anti-Choline Acetyltransferase/Chat Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P32738

### Technical Details

Immunogen	E.coli-derived rat Choline Acetyltransferase/Chat recombinant protein (Position: E19-D612).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm

	Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 7 Publications Citing This Product

1. PubMed ID: 10.1186/s12872-016-0375-3, Autonomic remodeling may be responsible for decreased incidence of aortic dissection in STZ-induced diabetic rats via down-regulation of matrix metalloprotease 2
2. PubMed ID: 10.1111/j.1527-3458.2007.00031.x, Protection of Mouse Brain from Aluminum<sup>3+</sup>induced Damage by Caffeic Acid
3. PubMed ID: 10.1111/cns.12187, Activation of Muscarinic Receptors Protects against Retinal Neurons Damage and Optic Nerve Degeneration In Vitro and In Vivo Models

Visit [bosterbio.com/anti-choline-acetyltransferase-chat-picoband-trade-antibody-a01192-5-boster.html](http://bosterbio.com/anti-choline-acetyltransferase-chat-picoband-trade-antibody-a01192-5-boster.html) to see all 7 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Choline Acetyltransferase/Chat Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.