

Anti-Fibrinogen Beta / FGB Antibody

Catalog Number: A01204

About FGB

Seems to phosphorylate critical substrates that regulate the G1/S phase transition and/or DNA replication. Can phosphorylates MCM2 and MCM3.

Sato N., EMBO J. 16:4340-4351(1997).

Hess G.F., Gene 211:133-140(1998).

Jiang W., Proc. Natl. Acad. Sci. U.S.A. 94:14320-14325(1997).

Overview

Product Name	Anti-Fibrinogen Beta / FGB Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Fibrinogen Beta / FGB Antibody catalog # A01204. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02675

Technical Details

Immunogen	Synthesized peptide derived from human MAP2 (Phospho-Thr1616)
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

Suggested Dilutions

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

WB: 1:500-1:2000

ICC/IF: 1:50-1:200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibrinogen Beta / FGB Antibody