

Anti-Vinculin/VCL Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A01207-1-Fluoro550

About VCL

In mammalian cells, vinculin is a membrane-cytoskeletal protein in focal adhesion plaques that is involved in linkage of integrin adhesion molecules to the actin cytoskeleton. It is mapped to 10q22.2. Vinculin is a cytoskeletal protein associated with cell-cell and cell-matrix junctions, where it is thought to function as one of several interacting proteins involved in anchoring F-actin to the membrane. Defects in VCL are the cause of cardiomyopathy dilated type 1W. Dilated cardiomyopathy is a disorder characterized by ventricular dilation and impaired systolic function, resulting in congestive heart failure and arrhythmia. Multiple alternatively spliced transcript variants have been found for this gene, but the biological validity of some variants has not been determined.

Overview

Product Name	Anti-Vinculin/VCL Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P18206

Technical Details

Immunogen	E.coli-derived human Vinculin/VCL recombinant protein (Position: M1-E149).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

3 Publications Citing This Product

1. PubMed ID: 10.1016/j.cbi.2016.03.013, Toxicity analysis of various Pluronic F-68-coated carbon nanotubes on mesenchymal stem cells
2. PubMed ID: 10.1007/s12010-015-1862-1, Monitoring of TGF-beta 1-Induced Human Lung Adenocarcinoma A549 Cells Epithelial-Mesenchymal Transformation Process by Measuring Cell Adhesion Force with a Microfluidic Device
3. PubMed ID: 32247806, Dong Y,Chen Y,Zhu D,Shi K,Ma C,Zhang W,Rocchi P,Jiang L,Liu X.Self-assembly of amphiphilic phospholipid peptide dendrimer-based nanovectors for effective delivery of siRNA therapeutics in prostate cancer therapy.J Control Release.2020 Jun 10;322:416-425.d

Visit bosterbio.com/anti-vinculin-vcl-picoband-trade-antibody-a01207-1-boster.html to see all 3 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Vinculin/VCL Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.