

Anti-P2X purinoceptor 7 P2RX7 Antibody

Catalog Number: A01208-1


Overview

Product Name	Anti-P2X purinoceptor 7 P2RX7 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-P2X purinoceptor 7 P2RX7 Antibody catalog # A01208-1. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99572

Technical Details

Immunogen	Recombinant fusion protein of human P2RX7(NP_002553.3).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000

Anti-P2X purinoceptor 7 P2RX7 Antibody (A01208-1) Images



Western blot analysis of extracts of Rat lung, using P2RX7 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 90s.

1 Publications Citing This Product

1. PubMed ID: 32783373, Zou Z,Wang Q,Zhou M,Li W,Zheng Y,Li F,Zheng S,He Z.Protective effects of P2X7R antagonist in sepsis-induced acute lung injury in mice via regulation of circ_0001679 and circ_0001212 and downstream Pln, Cdh2, and Npr13 expression.J Gene Med.2020 Dec;22(12)

Visit bosterbio.com/anti-p2x-purinoceptor-7-p2rx7-antibody-a01208-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-P2X purinoceptor 7 P2RX7 Antibody

For Research Use Only. Not for use in diagnostic procedures.