

## Anti-TGFBI Antibody

Catalog Number: A01218


### Overview

Product Name	Anti-TGFBI Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TGFBI Antibody catalog # A01218. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15582


### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200

## Anti-TGFBI Antibody (A01218) Images



Western blot analysis of TGFBI on mouse colon (1) and mouse eye (2) tissue lysate using anti-TGFBI antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue using anti-TGFBI antibody. Counter stained with hematoxylin.

## 1 Publications Citing This Product

1. PubMed ID: 32319604, Jia X, Wang H, Li Z, Yan J, Guo Y, Zhao W, Gao L, Wang B, Jia Y. HER4 promotes the progression of colorectal cancer by promoting epithelial-mesenchymal transition. Mol Med Rep. 2020 Apr;21(4):1779-1788. doi:10.3892/mmr.2020.10974. Epub 2020 Feb 4. PMID:32319604; PMCID

Visit [bosterbio.com/anti-tgfb1-antibody-a01218-boster.html](https://bosterbio.com/anti-tgfb1-antibody-a01218-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TGFBI Antibody

For Research Use Only. Not for use in diagnostic procedures.