

Anti-ABCG1 Antibody

Catalog Number: A01227


Overview

Product Name	Anti-ABCG1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ABCG1 Antibody catalog # A01227. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P45844


Technical Details

Immunogen	Synthesized peptide derived from ABCG1 . at AA range: 560-640
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-300 ELISA 1:10000 IF 1:50-200


Anti-ABCG1 Antibody (A01227) Images



Western blotting validation for Anti-ABCG1 Antibody A01227
Western Blot (WB) analysis of extracts from K562 cells



Immunohistochemistry validation of ABCG1 using Anti-ABCG1 Antibody (A01227). Immunohistochemistry (IHC) analysis of paraffin-embedded Rat Brain



Immunohistochemistry validation of ABCG1 using Anti-ABCG1 Antibody (A01227). Immunohistochemistry (IHC) analysis of paraffin-embedded Rat Brain

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ABCG1 Antibody

For Research Use Only. Not for use in diagnostic procedures.