

## Anti-NFkB-p100 (S869) NFKB2 Antibody

Catalog Number: A01228S869

### Overview


Product Name	Anti-NFkB-p100 (S869) NFKB2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NFkB-p100 (S869) NFKB2 Antibody catalog # A01228S869. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00653

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 840-890 of Human NFkB-p100.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-NFkB-p100 (S869) NFKB2 Antibody (A01228S869) Images

---



Immunohistochemistry (IHC) analyzes of NFkB-p100 (S869) pAb in paraffin-embedded human breast carcinoma tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NFkB-p100 (S869) NFKB2 Antibody

For Research Use Only. Not for use in diagnostic procedures.