

Anti-Tsg 101 (E303) Antibody

Catalog Number: A01233


Overview

Product Name	Anti-Tsg 101 (E303) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tsg 101 (E303) Antibody catalog # A01233. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99816


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 268-320 of Human Tsg 101.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200


Anti-Tsg 101 (E303) Antibody (A01233) Images



Western blot (WB) analysis of Tsg 101 (E303) polyclonal antibody in extracts from HEK293A cells.



Immunohistochemistry (IHC) analyzes of Tsg 101 (E303) pAb in paraffin-embedded human breast carcinoma tissue at 1:50, showing nucleus/cell membrane and cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.



Immunohistochemistry (IHC) analyzes of Tsg 101 (E303) pAb in paraffin-embedded human kidney carcinoma tissue at 1:50, showing cytoplasmic and nucleus staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

1 Publications Citing This Product

1. PubMed ID: -, Pei-pei Fang, Chen-wei Pan, Wei Lin, Jie Li, Shan-shan Huang, Guang-yao Zhou, Wen-jun Du, Qiang Li, "ASK1 Enhances Angiotensin II-Induced Liver Fibrosis In Vitro by Mediating Endoplasmic Reticulum Stress-Dependent Exosomes", Mediators of Inflammation, vol.2020, Art

Visit bosterbio.com/anti-tsg-101-e303-antibody-a01233-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tsg 101 (E303) Antibody

For Research Use Only. Not for use in diagnostic procedures.