

## Anti-Dysferlin Antibody

Catalog Number: A01234-2


### Overview

Product Name	Anti-Dysferlin Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Dysferlin Antibody catalog # A01234-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human,Mouse.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75923


### Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200

## Anti-Dysferlin Antibody (A01234-2) Images



Western blot analysis of Dysferlin on human skeletal muscle lysates using anti-Dysferlin antibody at 1/1,000 dilution.



ICC staining Dysferlin in SW480 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dysferlin Antibody

For Research Use Only. Not for use in diagnostic procedures.