

## Anti-USP7 Antibody

Catalog Number: A01239


### Overview

Product Name	Anti-USP7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-USP7 Antibody catalog # A01239. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q93009


### Technical Details

Immunogen	Recombinant protein within human USP7 aa 100-200.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50 FC: 1:50-1:100

## Anti-USP7 Antibody (A01239) Images



Western blot analysis of USP7 on 293T cell using anti-USP7 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-USP7 antibody. Counter stained with hematoxylin.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-USP7 Antibody

For Research Use Only. Not for use in diagnostic procedures.