

Anti-Apolipoprotein A-V APOA5 Antibody

Catalog Number: A01242-1


Overview

Product Name	Anti-Apolipoprotein A-V APOA5 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Apolipoprotein A-V APOA5 Antibody catalog # A01242-1. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6Q788

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 11-55 of Human APOA5
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Suggested Dilutions	WB: 1:500-1:1000

Anti-Apolipoprotein A-V APOA5 Antibody (A01242-1) Images



Western blot (WB) analysis of APOA5 polyclonal antibody at 1:500 dilution Lane1:HCT116 whole cell lysate(40ug) Lane2:The Brain tissue lysate of Mouse(40ug) Lane3:The Brain tissue lysate of Rat(40ug) Lane4:C6 whole cell lysate(40ug) Lane5:A549 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein A-V APOA5 Antibody

For Research Use Only. Not for use in diagnostic procedures.