

Anti-IRAK4 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A01247-3-Fluoro488

About IRAK4

IRAK-4(interleukin-1 receptor-associated kinase 4), also called REN64, in the IRAK family, is a protein kinase involved in signaling innate immune responses from Toll-like receptors. It also supports signaling from T-cell receptors. Scott(2002) mapped the REN64/IRAK4 gene to chromosome 12 based on similarity between the REN64 sequence(GenBank AF155118) and a chromosome 12 clone(GenBank AC093012). Functional analysis by Li et al.(2002) determined that IRAK4, like IRAK1 and Pelle, has auto- and cross-phosphorylation kinase activity. Precipitation and binding analyses showed weak interaction between IRAK4 and IRAK1, but IRAK4 did not interact with other IRAK family members. Overexpressed IRAK4 interacted with MYD88 and TRAF6 and activated mitogen-activated protein kinase(MAPK) and nuclear factor kappa-B(NFKB) pathways.

Overview

Product Name	Anti-IRAK4 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NWZ3

Technical Details

Immunogen	E.coli-derived human IRAK4 recombinant protein (Position: N419-S460).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IRAK4 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.