

Anti-HIF-2 Alpha EPAS1 Antibody

Catalog Number: A01248


Overview

Product Name	Anti-HIF-2 Alpha EPAS1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HIF-2 Alpha EPAS1 Antibody catalog # A01248. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99814


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-HIF-2 Alpha EPAS1 Antibody (A01248) Images



Western blot analysis of HIF-2 alpha on SiHa cell lysate using anti-HIF-2 alpha antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded rat lung tissue using anti-HIF-2 alpha antibody. Counter stained with hematoxylin.

1 Publications Citing This Product

1. PubMed ID: 33155655, Wang P, Meng QQ, Wang W, Zhang S, Xiong X, Qin S, Zhang J, Li A, Liu Z. Icaritin inhibits the inflammation through down-regulating NF-kappaB/HIF-2alpha signal pathways in chondrocytes. Biosci Rep. 2020 Nov 6; BSR20203107. doi:10.1042/BSR 20203 107. Epub ahead of print. PMID:3

Visit bosterbio.com/anti-hif-2-alpha-epas1-antibody-a01248-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HIF-2 Alpha EPAS1 Antibody

For Research Use Only. Not for use in diagnostic procedures.