

## Anti-HIF-2-alpha/EPAS1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A01248-1-Fluoro488

### About EPAS1

HIF-2 alpha is also designated EPAS1 whose gene is mapped to 2p21-p16. The predicted mouse protein is 88% identical to human EPAS1. The human EPAS1 gene contains 15 exons and spans at least 120 kb. The positions of the introns within the genomic region encoding the N-terminal bHLH-PAS domains of EPAS1 and AHR are similar, suggesting that the 5-prime ends of the 2 genes may have arisen from a gene duplication event. Moreover, the predicted protein shares 48% sequence identity with HIF1-alpha, a bHLH-PAS transcription factor that induces EPO gene expression in cultured cells in response to hypoxia. Like HIF1A, EPAS1 binds to and activates transcription from the HIF1A response element derived from the 3-prime flanking region of the EPO gene. EPAS1 is predominantly expressed in highly vascularized tissues of adult humans and in endothelial cells of the mouse adult and embryo. Furthermore, EPAS1 may represent an important regulator of vascularization, perhaps involving the regulation of endothelial cell gene expression in response to hypoxia. HIF2A is expressed at relatively higher levels in villus sections of placenta and in lung samples compared with other tissues examined. In addition, The variation in EPAS1 influences the relative contribution of aerobic and anaerobic metabolism and hence the maximum sustainable metabolic power for a given event duration.

### Overview

Product Name	Anti-HIF-2-alpha/EPAS1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q99814

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human HIF-2-alpha/EPAS1, which shares 94.4% amino acid (aa) sequence identity with both mouse and rat HIF-2-alpha/EPAS1.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL

Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HIF-2-alpha/EPAS1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.