

Anti-Chromogranin A CHGA Antibody

Catalog Number: A01256-1


Overview

Product Name	Anti-Chromogranin A CHGA Antibody
Reactive Species	Human
Description	Boster Bio Anti-Chromogranin A CHGA Antibody catalog # A01256-1. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10645


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-Chromogranin A CHGA Antibody (A01256-1) Images



Western blot analysis of Chromogranin A on PC-12 cells lysates using anti-Chromogranin A antibody at 1/1,000 dilution.



ICC staining Chromogranin A in PC-3M cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 31732520, Glumac PM, Gallant JP, Shapovalova M, Li Y, Murugan P, Gupta S, Coleman IM, Nelson PS, Dehm SM, LeBeau AM. Exploitation of CD133 for the Targeted Imaging of Lethal Prostate Cancer. Clin Cancer Res. 2020 Mar 1;26(5):1054-1064. doi:10.1158/1078-0432.CCR-19-1659. Epub 20

Visit bosterbio.com/anti-chromogranin-a-chga-antibody-a01256-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Chromogranin A CHGA Antibody

For Research Use Only. Not for use in diagnostic procedures.